

COVER FOR JAR 400 ML

Réf. BCASPI85



COVER FOR JAR 400 ML

CA-MI

Réf. BCASPI85

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169347
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		0,200	
		Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min		-20	
		Temperature max	45
Humidity min		0	
		Humidity max	85



30/03/2022 23:33 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI85

All characteristics are given for information only and may be changed

2/2

HOLTEX

30 rue Jean de Guirmand
13290 Aix-en-Provence - France

www.gsh-med.fr

Contact :

+33(0)4 42 90 31 31 - info@gsh-med.fr