

"O" RING FOR JAR 400 ML

Réf. BCASPI88



"O" RING FOR JAR 400 ML**CA-MI****Réf. BCASPI88**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169378
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:34 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI88

All characteristics are given for information only and may be changed

2 / 2