

## KOCHER PROSTATIC RETRACTOR, 24CM

Réf. IE22300



## KOCHER PROSTATIC RETRACTOR, 24CM



Réf. IE22300

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016986
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,5	Width packaging (cm):	
Length (cm) :	25,0	Length packaging (cm) :	
Height (cm) :	6,7	Height packaging (cm) :	
Weight (kg) :	0,216	Weight packaging (kg) :	0,216
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE22300

All characteristics are given for information only and may be changed