

MINI-LANGENBECK RETRACTOR 16 CM - 20 X 6 MM

Réf. IE25216



MINI-LANGENBECK RETRACTOR 16 CM - 20 X 6 MM



Réf. **IE25216**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040017020
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,0	Width packaging (cm):	
Length (cm) :	16,4	Length packaging (cm) :	
Height (cm) :	2,2	Height packaging (cm) :	
Weight (kg) :	0,063	Weight packaging (kg) :	0,063
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference IE25216

All characteristics are given for information only and may be changed