

LANGENBECK RETRACTOR 22 CM, 30 X 7 MM

Réf. IE25312



LANGENBECK RETRACTOR 22 CM, 30 X 7 MM



Réf. IE25312

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040084916
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	2,0	Width packaging (cm):	
Length (cm) :	22,5	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,074	Weight packaging (kg) :	0,074
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:36 - This article is available on our website : www.gsh-med.fr with the following reference IE25312

All characteristics are given for information only and may be changed