

STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK



Réf. OSA003_N



STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK



Réf. OSA003_N

- Rapid assembly system
- Square stainless steel tube frame
- Adjustable headboard with automatic high resistance rack
- Maximum weight: 150 kg
- Headboard and seat covered with 5 cm thick very high density foam (30 kg/m³), covered by fire-resistant washable plasticised fabric (fire resistance class: M1)
- Removable roll holder (for 50 cm wide rolls) that can be located at head or feet

OPTIONS

Réf. OTA0047 Heel support without clamps (pair)

Réf. OTA0041 Leg holders without clamps (pair)

Réf. OTA0044 Fixing clamps (pair)

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040165844
Guarantee (month):	36	Customs Code	9402900000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	63,0	Width packaging (cm):	15,0
Length (cm) :	190,0	Length packaging (cm) :	190,0
Height (cm) :	70,0	Height packaging (cm) :	68,0
Weight (kg) :	24,500	Weight packaging (kg) :	27,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:38 - This article is available on our website : www.gsh-med.fr with the following reference OSA003_N

All characteristics are given for information only and may be changed