

## EXAMINATION COUCH BLACK - EPOXY COATED STEEL (QUICK ASSEMBLING), BLACK

Réf. OSA003E\_N



# EXAMINATION COUCH BLACK - EPOXY COATED STEEL(QUIK ASSEMBLING), BLACK



Réf. OSA003E\_N

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040168135
Guarantee (month):	36	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	63,0	Width packaging (cm):	16,0
Length (cm) :	190,0	Length packaging (cm) :	190,0
Height (cm) :	70,0	Height packaging (cm) :	67,0
Weight (kg) :	25,000	Weight packaging (kg) :	28,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



05/09/2022 21:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OSA003E\_N