

## STAINLESS STEEL SCREEN 2 PANELS WITH STRAINED WHITE CURTAINS



Réf. OS0062



# STAINLESS STEEL SCREEN 2 PANELS WITH STRAINED WHITE CURTAINS



Réf. **OS0062**

- Panels (155 x 61 cm) with fixed castors without brakes
- Stainless steel tube frame
- Washable plastic curtains
- Fire resistance class: M3
- 4 double castors. Diameter: 60 mm

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040162621
Guarantee (month):	36	Customs Code	9403899000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	64,0
Length (cm) :		Length packaging (cm) :	177,0
Height (cm) :		Height packaging (cm) :	9,0
Weight (kg) :	6,000	Weight packaging (kg) :	7,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



18/10/2022 21:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OS0062

All characteristics are given for information only and may be changed