

FRAZIER CANNULA 18 CM X 4 MM

Réf. IC41104



FRAZIER CANNULA 18 CM X 4 MM



Réf. IC41104

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040011516
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,8	Width packaging (cm):	
Length (cm) :	20,5	Length packaging (cm) :	
Height (cm) :	1,8	Height packaging (cm) :	
Weight (kg) :	0,032	Weight packaging (kg) :	0,032
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:12 - This article is available on our website : www.gsh-med.fr with the following reference IC41104

All characteristics are given for information only and may be changed