

## SEMI-AUTOMATIC DEFIBRILLATOR - PRO-200J - SAVER ONE PRO - LCD SCREEN

Réf. SJSAVE10

Saver One®



GARANTIE  
**6**  
ANS  
YEARS  
GUARANTEE

CE  
0051

DISPOSITIF MEDICAL CLASSE  
MEDICAL DEVICE IIb

## SEMI-AUTOMATIC DEFIBRILLATOR - PRO-200J - SAVER ONE PRO - LCD SCREEN



Réf. SJSAVE10

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0051	Barcode	3661040128726
Guarantee (month):	72	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	2,000	Weight packaging (kg) :	4,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:13 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SJSAVE10

All characteristics are given for information only and may be changed