

EUROPE TO UK POWER ADAPTER

Réf. OTA0909



EUROPE TO UK POWER ADAPTER

Réf. **OTA0909**

FEATURES			
Medical device	No		
Marquage	Barcode	3661040039893	
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	10,0	Width packaging (cm):	
Length (cm) :	14,0	Length packaging (cm) :	
Height (cm) :	5,0	Height packaging (cm) :	
Weight (kg) :	0,045	Weight packaging (kg) :	0,045
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:12 - This article is available on our website : www.gsh-med.fr with the following reference OTA0909

All characteristics are given for information only and may be changed