

## EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING), BLUE LAVANDER



Réf. OSA003\_L



# EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING),BLUE LAVANDER



Réf. OSA003\_L

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040196152
Guarantee (month):	36	Customs Code	9402900000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	63,0	Width packaging (cm):	15,0
Length (cm) :	190,0	Length packaging (cm) :	190,0
Height (cm) :	70,0	Height packaging (cm) :	68,0
Weight (kg) :	24,500	Weight packaging (kg) :	27,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/10/2022 12:38 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OSA003\_L