

ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE

Réf. SBRAN08_O



ALUMINIUM FOLDABLE STRETCHER IN 2 (LENGHT AND LARG) + POUCH, ORANGE



Réf. **SBRAN08_O**

Folding stretcher with balancing system

- High resistance aluminium structure and waterproof fabric, fire resistant and anti-cracking
- This balancing stretcher folds in length and width
- Lightweight, easy to carry
- Dimensions (L x W x H): 208 X 55 x 13 cm, unfolded; 104 x 17 x 9 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Single colour: orange
- Delivered with a stretcher bag

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094632
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	54,5	Width packaging (cm):	20,0
Length (cm) :	208,0	Length packaging (cm) :	10,0
Height (cm) :	13,5	Height packaging (cm) :	10,0
Weight (kg) :	4,700	Weight packaging (kg) :	5,600
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:20 - This article is available on our website : www.gsh-med.fr with the following reference SBRAN08_O

All characteristics are given for information only and may be changed

2/2

HOLTEX

30 rue Jean de Guiramand
13290 Aix-en-Provence - France

www.gsh-med.fr

Contact :

+33(0)4 42 90 31 31 - info@gsh-med.fr