

## ALUMINIUM ALLOY FOLDAWAY STRETCHER, ORANGE

Réf. SBRAN09\_O



## ALUMINIUM ALLOY FOLDAWAY STRETCHER, ORANGE



Réf. **SBRAN09\_0**

- Folding stretcher with balancing system
- High resistance aluminium structure and waterproof fabric, fire resistant and anti-cracking
  - This balancing stretcher folds in length
  - Lightweight, easy to carry
  - Dimensions (L x W x H): 215 X 54.5 x 13 cm, unfolded; 182 x 17 x 9 cm, folded
  - Weight: 8 kg
  - Maximum load carried: 159 kg
  - Single colour: orange

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094649
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	54,5	Width packaging (cm):	20,0
Length (cm) :	207,0	Length packaging (cm) :	211,0
Height (cm) :	13,0	Height packaging (cm) :	10,0
Weight (kg) :	5,000	Weight packaging (kg) :	6,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:20 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SBRAN09\_0