

HE.CO STOP COMPRESSION CUSHION BANDAGE

Réf. SPSECU86



HE.CO STOP COMPRESSION CUSHION BANDAGE



Réf. **SPSECU86**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040116945
Guarantee (month):	0	Customs Code	30059100
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	11,2	Width packaging (cm):	16,0
Length (cm) :	17,0	Length packaging (cm) :	19,0
Height (cm) :	4,3	Height packaging (cm) :	4,0
Weight (kg) :	0,070	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:14 - This article is available on our website : www.gsh-med.fr with the following reference SPSECU86

All characteristics are given for information only and may be changed