

BECKMANN RETRACTOR (4X4 TEETH) 13 CM

Réf. IE13413



BECKMANN RETRACTOR (4X4 TEETH) 13 CM



Réf. IE13413

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016610
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,2	Width packaging (cm):	
Length (cm) :	14,5	Length packaging (cm) :	
Height (cm) :	1,8	Height packaging (cm) :	
Weight (kg) :	0,064	Weight packaging (kg) :	0,064
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:41 - This article is available on our website : www.gsh-med.fr with the following reference IE13413

All characteristics are given for information only and may be changed