

BECKMANN RETRACTOR (4X4 DIENTES) 16 CM

Réf. IE13416



BECKMANN RETRACTOR (4X4 DIENTES) 16 CM



Réf. IE13416

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016627
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,3	Width packaging (cm):	
Length (cm) :	16,4	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,088	Weight packaging (kg) :	0,088
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference IE13416

All characteristics are given for information only and may be changed