

BECKMANN RETRACTOR (4X4 TEETH) 30 CM

Réf. IE13430



BECKMANN RETRACTOR (4X4 TEETH) 30 CM



Réf. IE13430

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040016641
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,7	Width packaging (cm):	
Length (cm) :	30,5	Length packaging (cm) :	
Height (cm) :	2,7	Height packaging (cm) :	
Weight (kg) :	0,220	Weight packaging (kg) :	0,220
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference IE13430

All characteristics are given for information only and may be changed