

STAINLESS STEEL SCREEN 4 PANELS WITH STRAINED WHITE CURTAINS- HOLTEX



Réf. OS0064



STAINLESS STEEL SCREEN 4 PANELS WITH STRAINED WHITE CURTAINS- HOLTEX



Réf. OS0064

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040168111
Guarantee (month):	36	Customs Code	9403899000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	64,0
Length (cm) :		Length packaging (cm) :	177,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :	10,600	Weight packaging (kg) :	13,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



18/10/2022 21:11 - This article is available on our website : www.gsh-med.fr with the following reference OS0064

All characteristics are given for information only and may be changed