

STAINLESS STEEL STIRRUPS FOR VOG ELECTRIC EXAMINATION COUCH

Réf. OV0013



STAINLESS STEEL STIRRUPS FOR VOG ELECTRIC EXAMINATION COUCH



Réf. OV0013

FEATURES			
Medical device	No		
Marquage		Barcode	3661040094342
Guarantee (month):	36	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:12 - This article is available on our website : www.gsh-med.fr with the following reference OV0013

All characteristics are given for information only and may be changed