

ARNEODO HOOK FOR PHACO EMULSIFICATION

Réf. IC60480



ARNEODO HOOK FOR PHACO EMULSIFICATION



Réf. IC60480

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040084046
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,5	Width packaging (cm):	1,5
Length (cm) :	12,1	Length packaging (cm) :	12,1
Height (cm) :	0,5	Height packaging (cm) :	0,5
Weight (kg) :	0,013	Weight packaging (kg) :	0,016
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:22 - This article is available on our website : www.gsh-med.fr with the following reference IC60480

All characteristics are given for information only and may be changed