

WEBER IRIS KNIFE STRAIGHT, 13CM

Réf. IC60500



WEBER IRIS KNIFE STRAIGHT, 13CM

Réf. IC60500

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040012766
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,6	Width packaging (cm):	
Length (cm) :	12,9	Length packaging (cm) :	
Height (cm) :	0,4	Height packaging (cm) :	
Weight (kg) :	0,050	Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference IC60500

All characteristics are given for information only and may be changed