

RAY PITUITARY CURETTE HORIZONTAL 90°, 4 MM.

Réf. IC70100



RAY PITUITARY CURETTE HORIZONTAL 90°, 4 MM.



Réf. **IC70100**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040012988
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,7	Width packaging (cm):	
Length (cm) :	18,0	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,030	Weight packaging (kg) :	0,030
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:22 - This article is available on our website : www.gsh-med.fr with the following reference IC70100

All characteristics are given for information only and may be changed