

MAURER RIB SHEAR RIGHT 36 CM

Réf. IC32836



MAURER RIB SHEAR RIGHT 36 CM



Réf. IC32836

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040011066
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	13,5	Width packaging (cm):	
Length (cm) :	36,5	Length packaging (cm) :	
Height (cm) :	1,2	Height packaging (cm) :	
Weight (kg) :	0,463	Weight packaging (kg) :	0,463
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:08 - This article is available on our website : www.gsh-med.fr with the following reference IC32836

All characteristics are given for information only and may be changed