

Réf. DNACCE21



CE
0123



Réf. DNACCE21

- Flexible, several sizes to suit all morphologies
- Variety of cable lengths
- Reusable sensors

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0123	Barcode	3661040124315
Guarantee (month):	3	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	300,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:00 - This article is available on our website : www.gsh-med.fr with the following reference DNACCE21