

## SLIDING SOCKETS FOR VOG ELECTRIC EXAMINATION COUCH

Réf. OV0015



# SLIDING SOCKETS FOR VOG ELECTRIC EXAMINATION COUCH



Réf. OV0015

FEATURES			
Medical device	No		
Marquage		Barcode	3661040094366
Guarantee (month):	36	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,350
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:12 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference OV0015

All characteristics are given for information only and may be changed