

## PEN LIGHT

Réf. DLE1020



## PEN LIGHT



Réf. DLE1020

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3701126702935
Guarantee (month):	12	Customs Code	85131010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,8	Width packaging (cm):	
Length (cm) :	16,5	Length packaging (cm) :	
Height (cm) :	1,8	Height packaging (cm) :	
Weight (kg) :	0,039	Weight packaging (kg) :	0,039
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:53 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DLE1020