

COMBAT APPLICATION TOURNIQUET, BLACK COLOR

Réf. STGARR01



COMBAT APPLICATION TOURNIQUET, BLACK COLOR



Réf. STGARR01

| FEATURES | | | |
|----------------------------|-----------|-------------------------|------------|
| Medical device | Yes | Classe I | |
| Marquage | CE | Barcode | |
| Guarantee (month): | 12 | Customs Code | 3006100000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | 0,078 | Weight packaging (kg) : | 0,078 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



14/02/2023 12:57 - This article is available on our website : www.gsh-med.fr with the following reference STGARR01

All characteristics are given for information only and may be changed