

BOWMAN SPECULUM 7 CM

Réf. IB10007



BOWMAN SPECULUM 7 CM



Réf. **IB10007**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040008400
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,8	Width packaging (cm):	6,8
Length (cm) :	20,0	Length packaging (cm) :	20,0
Height (cm) :	2,0	Height packaging (cm) :	2,0
Weight (kg) :	0,005	Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:54 - This article is available on our website : www.gsh-med.fr with the following reference IB10007

All characteristics are given for information only and may be changed