

PLEY SPECULUM 7 CM

Réf. IB10107



PLEY SPECULUM 7 CM



Réf. IB10107

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040008417
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,9	Width packaging (cm):	
Length (cm) :	7,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,017	Weight packaging (kg) :	0,017
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



22/08/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference IB10107

All characteristics are given for information only and may be changed