

GUYON CATHETERS WITH GUIDE N°55

Réf. IB03255



GUYON CATHETERS WITH GUIDE N°55



Réf. **IB03255**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040036212
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:49 - This article is available on our website : www.gsh-med.fr with the following reference IB03255

All characteristics are given for information only and may be changed