

## WHITE PLASTIC HARD EARTIPS IDEAL/ZAPHYR/CONSULTO (PAIR)

Réf. DS80015



# WHITE PLASTIC HARD EARTIPS IDEAL/ZAPHYR/CONSULTO (PAIR)



Réf. **DS80015**

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | No        |                         |               |
| Marquage                   |           | Barcode                 | 3661040001685 |
| Guarantee (month):         | 12        | Customs Code            |               |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 1,0       | Width packaging (cm):   |               |
| Length (cm) :              | 1,0       | Length packaging (cm) : |               |
| Height (cm) :              | 1,0       | Height packaging (cm) : |               |
| Weight (kg) :              | 0,091     | Weight packaging (kg) : | 0,091         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



30/03/2022 14:58 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference DS80015

All characteristics are given for information only and may be changed