

MARTIN FLEXIBLE UTERINE SOUND 30 CM

Réf. IB08603



MARTIN FLEXIBLE UTERINE SOUND 30 CM



Réf. IB08603

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040061344
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,5	Width packaging (cm):	
Length (cm) :	38,0	Length packaging (cm) :	
Height (cm) :	0,2	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,070
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:50 - This article is available on our website : www.gsh-med.fr with the following reference IB08603

All characteristics are given for information only and may be changed