

CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E

Réf. DPOXY04



CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E



Réf. **DPOXY04**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040163833
Guarantee (month):	6	Customs Code	9018909911
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :	0,059	Weight packaging (kg) :	0,001
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:12 - This article is available on our website : www.gsh-med.fr with the following reference DPOXY04

All characteristics are given for information only and may be changed