

## RECHARGEABLE BATTERY PACK, INTERNAL FOR KERN FOB(215X215 MM)&KERN MCB

***KERN***

Réf. EKFOB01



# RECHARGEABLE BATTERY PACK, INTERNAL FOR KERN FOB(215X215 MM)&KERN MCB



Réf. **EKFOB01**

| FEATURES                   |       |                         |               |
|----------------------------|-------|-------------------------|---------------|
| Medical device             | No    |                         |               |
| Marquage                   |       | Barcode                 | 4045761067134 |
| Guarantee (month):         | 0     | Customs Code            |               |
| Sterile                    | No    |                         |               |
| INFORMATIONS               |       | Main unit               |               |
|                            |       | Packaged unit           |               |
| Width (cm) :               |       | Width packaging (cm):   | 8,5           |
| Length (cm) :              |       | Length packaging (cm) : | 11,5          |
| Height (cm) :              |       | Height packaging (cm) : | 3,5           |
| Weight (kg) :              |       | Weight packaging (kg) : | 0,200         |
| SPECIFIC STORAGE CONDITION |       |                         |               |
| Temperature min            | (N/A) | Temperature max         | (N/A)         |
| Humidity min               | (N/A) | Humidity max            | (N/A)         |



30/03/2022 15:03 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKFOB01

All characteristics are given for information only and may be changed