

ADENOID SET

Réf. IBOI330



ADENOID SET



Réf. IBOI330

FEATURES

Medical device Yes Classe I

Marquage CE

Barcode 3661040082455

Guarantee (month): 24

Customs Code 9018.9090

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) :

Width packaging (cm):

Length (cm) :

Length packaging (cm) :

Height (cm) :

Height packaging (cm) :

Weight (kg) :

Weight packaging (kg) : 2,147

SPECIFIC STORAGE CONDITION

Temperature min (N/A)

Temperature max (N/A)

Humidity min (N/A)

Humidity max (N/A)



30/03/2022 15:21 - This article is available on our website : www.gsh-med.fr with the following reference IBOI330

All characteristics are given for information only and may be changed

ADENOID SET

Réf. IBOI330

COMPOSITION

Code	Descriptive	Quantity
IB60419	Killian Tonsil Scalpel 19 cm	1
IC08518	Metzemaun Scissors Curved 18 cm	1
IC09511	Iris Scissors Curved 11.5 cm	1
IC35510	Beckmann Knife for Adenoid growths 22 cm x 10 mm	1
IC35512	Beckmann Knife for Adenoid growths 22 cm x 12 mm	1
IC35514	Beckmann Knife for Adenoid growths 22 cm x 14 mm	1
IC35516	Beckmann Knife for Adenoid growths 22 cm x 16 mm	1
IC46127	Yankauer Cannula 27 cm	1
IC70110	Blondiau Curette (E.N.T), 20mm x 24cm	2
ID03602	Doyen Mouth-Gag Adult, 13.5 cm	1
ID05403	Scalpel Handle N°3	1
IE32100	Hurd Dissector 22 cm	1
IF25110	Loops By Dozen 0.3 mm	1
AI28146	S/S Box 28 x 14 x 6 cm	1
IO08413	Jennings Mouth-Gag 13 cm	1
IP44240	Chatelier adenoid forceps 22 cm	1
IP46419	Colver Tonsil Forceps 19 cm	1
IP46819	Negus Forceps very Curved 19 cm	1
IP47403	Hartmann Punch Forceps 7 mm	1
IS18428	Bruening Tonsil Snare, 28cm	1