


RECHARGEABLE BATTERY FOR MCC/MPE/MPD/MPC/MBC

KERN

Réf. EK MBC08



RECHARGEABLE BATTERY FOR MCC/MPE/MPD/MPC/MBC



Réf. **EKMBC08**

FEATURES			
Medical device	No		
Marquage	AUCUN	Barcode	3661040173535
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,184
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:28 - This article is available on our website : www.gsh-med.fr with the following reference EKMBC08

All characteristics are given for information only and may be changed