

KERN MBC ELECTRONIC BABY SCALE - 15 KGS - MEDICAL CLASS III

Réf. EKMB01_CIII

KERN



CE
0297



KERN MBC ELECTRONIC BABY SCALE - 15 KGS - MEDICAL CLASS III



Réf. EKMBC01_CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297		Barcode 3661040054995
Guarantee (month):	24		Customs Code
Sterile	No		
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:03 - This article is available on our website : www.gsh-med.fr with the following reference EKMBC01_CIII

KERN MBC ELECTRONIC BABY SCALE - 15 KGS - MEDICAL CLASS III



Réf. EKMBC01_CIII

COMPOSITION

Code	Descriptive	Quantity
EKMBC01	Baby scale MBC	1
EKCERTI05	Homologation certificate 965-128_ (Class III Medical)	1