

KERN MBC ELECTRONIC BABY SCALE - 20 KGS - MEDICAL CLASS III

Réf. EKMBC02_CIII

KERN



CE
0297



KERN MBC ELECTRONIC BABY SCALE - 20 KGS - MEDICAL CLASS III



Réf. EKMBC02_CIII

FEATURES					
Medical device	Yes	Classe Im (fonction de mesurage)		Accuracy measurement	Classe III
Marquage	CE0297		Barcode	3661040055008	
Guarantee (month):	24		Customs Code		
Sterile	No				
INFORMATIONS	Main unit		Packaged unit		
Width (cm) :		Width packaging (cm):			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :			
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)		Temperature max	(N/A)
Humidity min		(N/A)		Humidity max	(N/A)



30/03/2022 15:03 - This article is available on our website : www.gsh-med.fr with the following reference EKMBC02_CIII

All characteristics are given for information only and may be changed

KERN MBC ELECTRONIC BABY SCALE - 20 KGS - MEDICAL CLASS III



Réf. EK MBC02_CIII

COMPOSITION

Code	Descriptive	Quantity
EKMBC02	Baby Scale Kern MBC, 20Kg	1
EKCERTI05	Homologation certificate 965-128_ (Class III Medical)	1