

KERN MBC ELECTRONIC BABY SCALE - 6-15 KG WITH HEIGHT ROD - CLASS III

Réf. EKMB06-CIII

KERN



CE
0297



KERN MBC ELECTRONIC BABY SCALE - 6-15 KG WITH HEIGHT ROD - CLASS III



Réf. EKMBC06-CIII

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe III
Marquage	CE0297	Barcode	3661040063195
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	89,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	5,500	Weight packaging (kg) :	8,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:04 - This article is available on our website : www.gsh-med.fr with the following reference EKMBC06-CIII

KERN MBC ELECTRONIC BABY SCALE - 6-15 KG WITH HEIGHT ROD - CLASS III



Réf. EKMBC06-CIII

COMPOSITION

Code	Descriptive	Quantity
EKMBC01	Baby scale MBC	1
EKMBC03	height meter for Baby Scale MBC	1
EKCERTI05	Homologation certificate 965-128_ (Class III Medical)	1