

CONGENITAL CATARACT SURGERY SET

Réf. IBOI413



CONGENITAL CATARACT SURGERY SET



Réf. IBOI413

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040082769
Guarantee (month):	24	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,383
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:24 - This article is available on our website : www.gsh-med.fr with the following reference IBOI413

All characteristics are given for information only and may be changed

CONGENITAL CATARACT SURGERY SET

Réf. IBOI413

COMPOSITION

Code	Descriptive	Quantity
IB10405	Barraquer-Colibri Speculum 5 cm	1
IC60200	Graefe Cataract Knife, 13cm	1
IC60201	Paufique Cataract Knife	2
IC92103	Rycroft Air and Liquid Injection Cannula 0.3 mm	1
IC92113	Bonnet Irrigation & Aspiration cannula angled 2 mm	1
IC92114	Barraquer Alpha-Chymotrypsin injection cannula bulbous type	1
AI18084	S/S Box 18 x 8 x 4 cm	1