

KERN ELECTRONIC HAND GRIP DYNAMOMETER MPA (10,20,40,80KG)

Réf. EKMAP02

KERN

GARANTIE
2
ANS
YEARS
GARANTIE



KERN ELECTRONIC HAND GRIP DYNAMOMETER MPA (10,20,40,80KG)



Réf. EKMAP02

FEATURES			
Medical device	No		
Marquage	Barcode	4045761095762	
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	26,5	Width packaging (cm):	
Length (cm) :	35,0	Length packaging (cm) :	
Height (cm) :	8,5	Height packaging (cm) :	
Weight (kg) :	0,300	Weight packaging (kg) :	1,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 15:03 - This article is available on our website : www.gsh-med.fr with the following reference EKMAP02

All characteristics are given for information only and may be changed