

ELITE BAGS ISOTHERMAL / DIABETES BAG -DIABETIC -BLUE JEAN (EB14.017)

 ELITE BAGS

Réf. MEDIAB03_B



ELITE BAGS ISOTHERMAL / DIABETES BAG -DIABETIC -BLUE JEAN (EB14.017)



Réf. MEDIAB03_B

FEATURES			
Medical device	No		
Marquage		Barcode	8436564542458
Guarantee (month):	24	Customs Code	4202.12.99.90
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,0	Width packaging (cm):	19,0
Length (cm) :	17,0	Length packaging (cm) :	10,5
Height (cm) :	10,0	Height packaging (cm) :	6,0
Weight (kg) :	0,012	Weight packaging (kg) :	0,248
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:02 - This article is available on our website : www.gsh-med.fr with the following reference MEDIAB03_B