

VAQUEZ-LAUBRY CLASSIC PULL ON CUFF WITH BLADDER/DARK BLUE/ADULT (M)

Réf. 330112



VAQUEZ-LAUBRY CLASSIC PULL ON CUFF WITH BLADDER/DARK BLUE/ADULT (M)



Réf. **330112**

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Straps fastening.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446031886
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,225
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:32 - This article is available on our website : www.gsh-med.fr with the following reference 330112

All characteristics are given for information only and may be changed