

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/ADULT (M)

Réf. 330012



VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/ADULT (M)



Réf. 330012

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446028541
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,190	Weight packaging (kg) :	0,190
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:24 - This article is available on our website : www.gsh-med.fr with the following reference 330012

All characteristics are given for information only and may be changed