

SPENGLER MAGISTER II STETHOSCOPE, DUAL, SILVER

Réf. 507251



SPENGLER MAGISTER II STETHOSCOPE, DUAL, SILVER



Réf. 507251

Stethoscope, binaural, complete: A mechanical listening device designed for listening to sounds from the heart and lungs from adult and paediatric patients.

VERY HIGH ACOUSTIC QUALITY

- Dual frequency floating membrane:
 - Light pressure to hear low frequencies (vascular and cardiac)
 - Firm pressure to hear high frequencies (pulmonary)
- High performance single tubing
- Aluminium binaural

COMFORT

- Non-chill ring
- Silicone earpieces
- Bell indicator
- New binaural design better adapted to human anatomy.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446029968
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	16,5
Length (cm) :		Length packaging (cm) :	32,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,430
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : www.gsh-med.fr with the following reference 507251

All characteristics are given for information only and may be changed