

IDEAL+ DOUBLE HEAD ADULT, BLUE

Réf. DH7073_B



IDEAL+ DOUBLE HEAD ADULT, BLUE



Réf. DH7073_B

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661041000076
Guarantee (month):	12	Customs Code	9018908400
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	25,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,160
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:29 - This article is available on our website : www.gsh-med.fr with the following reference DH7073_B