

IDEAL+ DOUBLE HEAD PEDIATRIC PINK

Réf. DH7073_RC



IDEAL+ DOUBLE HEAD PEDIATRIC PINK



Réf. DH7073_RC

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661041000106
Guarantee (month):	12	Customs Code	9018908400
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	12,0
Length (cm) :		Length packaging (cm) :	25,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:23 - This article is available on our website : www.gsh-med.fr with the following reference DH7073_RC

All characteristics are given for information only and may be changed

2/2