

IDEAL+ SINGLE HEAD ADULT RED

Réf. DH7081_R



IDEAL+ SINGLE HEAD ADULT RED



Réf. DH7081_R

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661041000137
Guarantee (month):	12	Customs Code	9018908400
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,150
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:23 - This article is available on our website : www.gsh-med.fr with the following reference DH7081_R

All characteristics are given for information only and may be changed

2/2