

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/RED/ADULT (M)

Réf. 330013



VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/RED/ADULT (M)



Réf. **330013**

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | | Barcode | 3700446008833 |
| Guarantee (month): | 6 | Customs Code | 9018901000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | | Weight packaging (kg) : | 0,190 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 02:22 - This article is available on our website : www.gsh-med.fr with the following reference 330013

All characteristics are given for information only and may be changed